

The ibidi product family is comprised of a variety of μ-Slides and μ-Dishes, which have all been designed for high-end microscopic analysis of fixed or living cells.

The glass bottom versions of the μ-Slides and μ-Dishes are especially designed for TIRF and single molecule applications.

The μ-Slide 2 well <sup>Ph+</sup> (Phase contrast plus) is an array of 2 square fields where cells can be cultivated and investigated with microscopical methods. The μ-Slide 2 well <sup>Ph+</sup> improves the optical quality of phase contrast microscopy. In contrast to the classic μ-Slide 2 well, the <sup>Ph+</sup> version provides a special plate in the center of the wells. This plate suppresses

the meniscus which is disturbing the phase contrast effect in normal open wells. Openings near the corners provide access to the wells for filling and aspirating liquids easily.

## Material

The glass bottom version of the μ-Slides are made of a standard μ-Slide but with a glass coverslip bottom. It is not possible to detach the bottom. The μ-Slides are not autoclavable since they are temperature stable only up to 80°C / 175°F.

### Optical Properties ibidi glass bottom

Refractive index $n_D$	1.523
Abbe number	55
Thickness	No. 1.5H (selected quality 170 μm, ± 5 μm)
Material	Schott borosilicate glass, D 263M

## Surface and coating

The μ-Slide 2 well <sup>Ph+</sup> glass bottom is manufactured with an uncoated glass coverslip. Washing steps (e.g. with PBS) before cell seeding can remove glass dust which is advantageous for direct cell growth on the surface.

Protein coatings increase direct cell growth of adherent cells. Specific coatings on glass are possible following this protocol:

- Prepare your coating solution according to the manufacturer's specifications or reference. Prepare your μ-Slide. Adjust the concentration to a coating area of 11.4 cm<sup>2</sup> and 1.5 ml.
- Apply 1.5 ml into the growth area. Make sure that the entire bottom is covered with liquid easily tilting or shaking the μ-Slide. Put on the lid and leave at room temperature for at least 30 minutes.
- Aspirate the solution and wash with the recommended protein dilution buffer. Optionally, let dry at room temperature. Attention, some coating proteins might degenerate when drying!

## Geometry

The μ-Slide 2 well <sup>Ph+</sup> glass bottom provides a standard slide format according to ISO 8037/1.

### Geometry of μ-Slide 2 well <sup>Ph+</sup> glass bottom

Number of wells	2
Dimensions of wells (w × l × h) in mm	21.2 × 23.3 × 3.0
Growth area per well	4.8 cm <sup>2</sup>
Coating area per well	11.4 cm <sup>2</sup>
Volume per well	1.5 ml
Liquid height	3.0 mm
Total height with lid	10.8 mm
Bottom matches coverslip	No. 1.5

## Seeding cells

- Trypsinize and count cells as usual. Dilute the cell suspension to the desired concentration. Depending on your cell type, application of a 5-11 × 10<sup>4</sup> cells/ml suspension should result in a confluent layer within 2-3 days.
- Apply 1.5 ml cell suspension into each well of the μ-Slide. Avoid shaking as this will result in inhomogeneous distribution of the cells.
- Cover the slide with the supplied lid. Incubate at 37°C and 5% CO<sub>2</sub> as usual.

Undemanding cells can be left in their seeding medium for up to three days and grow to confluence there. However, best results might be achieved when the medium is changed every 1-2 days. Carefully aspirate the old medium and replace it by 1.5 ml/well fresh medium.

**Tip:**

The day before seeding the cells we recommend placing the cell medium and the μ-Slide into the incubator for equilibration. This will prevent the liquid inside from emerging air bubbles over the incubation time.


**Solvents for Fixation and Staining**

Cells can be observed live or fixed directly in the μ-Slide preferably on an inverted microscope. The slide material is compatible to acids, alkalis, PFA, and silicone oil. Alcohols may be used for short term incubation (e.g. cell fixation). Acetone is not compatible. Further specifications can be found at [www.ibidi.com](http://www.ibidi.com).

For optimal results in fluorescence microscopy and storage of stained probes ibidi provides a mounting medium (50001) optimized for μ-Dishes and μ-Slides.

**Filling and Handling**




Fill the wells by using a standard pipet. Inject the cell suspension directly into one of the openings. Medium exchange is easily done by aspirating the entire volume and refilling using 1.5 ml per well.



**μ-Slide 2 well Ph+ selection guide**





**μ-Slide 2 well**

Standard open wells for maximum sample access. Meniscus disturbs the beam path. Good phase contrast quality only in the center of each well



**μ-Slide 2 well Ph+**

Special plate in the center of the wells suppresses meniscus formation. No meniscus – parallel beam path. For excellent phase contrast microscopy all over the wells.



**Immersion Oil**

When using oil immersion objectives, use only the immersion oils specified in the table. The use of a non-recommended oil could lead to the damage of the plastic material and the objective.

Company	Product	Ordering Number
ibidi	Immersion Oil	(ibidi) 50101
Zeiss	Immersol 518 F	(Zeiss) 444960
Zeiss	Immersol W 2010	(Zeiss) 444969
Leica	Immersion liquid	(Leica) 11513859

**μ-Slide 2 well Family**


The μ-Slide 2 well family is available as open well and as a Ph+ version. See table below for choosing your μ-Slide 2 well.

μ-Slide 2 well



Ordering Number	Treatment or Coating	Characteristics
80286	ibiTreat, sterile	hydrophilic, tissue culture treated
80282	Collagen IV, sterile	protein coating
80283	Fibronectin, sterile*	protein coating
80284	Poly-L-Lysine, sterile	biopolymer coating
80285	Poly-D-Lysine, sterile*	biopolymer coating
80281	uncoated, sterile	hydrophobic
80287	glass bottom	glass coverslip No. 1.5H (170 μm ±5 μm)

μ-Slide 2 well<sup>Ph+</sup>



Ordering Number	Treatment or Coating	Characteristics
80296	ibiTreat, sterile	hydrophilic, tissue culture treated
80292	Collagen IV, sterile	protein coating
80293	Fibronectin, sterile*	protein coating
80294	Poly-L-Lysine, sterile	biopolymer coating
80295	Poly-D-Lysine, sterile*	biopolymer coating
80291	uncoated, sterile	hydrophobic
80297	glass bottom	glass coverslip No. 1.5H (170 μm ±5 μm)

\* available on request only

**For research use only!**

Further technical specifications can be found at [www.ibidi.com](http://www.ibidi.com). For questions and suggestions please contact us by e-mail [info@ibidi.de](mailto:info@ibidi.de) or by telephone +49 (0)89/520 4617 0. All products are developed and produced in Germany.  
© ibidi GmbH, Am Klopferspitz 19, 82152 Martinsried, Germany.